




- Arch to operate, antagonist and centric
 - Material:
 - If teeth are stable: silicone (polyether or similar polysiloxanes)
 - If teeth are not stable: alginate (after removal or splinting)
 - Maximum extension, non-functionalised (anatomical)
- Reproduce the exact situation of the mouth at the time of surgery


- Extra-hard plaster model casting (no imperfections or bubbles)




- Fix the 3DMarkers to the bite according to the diagram above and glue them with a drop of cyanoacrylate
- Measure the bite on the model and if necessary shorten and fit it to the dimensions of the arch. It is advisable for the bite to reach the molars (not over)

- Arch to operate, antagonist and centric
 - Maximum extension, without holes
- Reproduce the exact situation of the mouth at the time of surgery

- High definition 3D print of the STL file




- Secure the Evobite with transparent and not radiopaque resin (orthodontic type) on the model, particularly compensating the space in line with the edentulous areas
- Keep the front 3DMarker more central and as near to the incisor edge as possible.



- Apply a layer of universal adhesive between the silicone and the bite and apply the silicone on the side of the bite that interlocks with the model
- Compensate with the silicone on the model corresponding to the edentulous area
- Fix the bite on the model taking care that you keep the 3DMarker more centred and as near to the incisor edge as possible


When hardened, position the bite on the patient and reline towards the antagonist with the same silicone supplied with the Evobite kit, taking care that you keep the bite steady (note: it is not a centric recording but a simple stabilisation key)




NOTE: THE COMPONENTS TO CONSTRUCT THE EVOBITE ARE DISPOSABLE



TRADITIONAL IMPRESSION



CASTING MODELS




VERTICAL DIMENSION AND OCCLUSAL RELATIONSHIP DETECTION AESTHETIC TEST



PRECISION DUPLICATION WITH TRANSPARENT RESIN




EVOBITE MOUNTING WITH SILICON ON THE MODEL WITH THE PROSTHESIS DUPLICATE MOUNTED



COMPATIBLE PROTOCOL

The RealGUIDE procedure is compatible with DUAL SCAN protocol (double scan) that consists of the insertion of radiopaque points in the prosthesis duplicate and double CBCT scan (patient with radiological template and radiological template only) to automatically overlap the scan of the prosthesis with the anatomy of the patient directly into the 3Diagnosis software. However, it is advisable to follow the standard procedure and use the Evobite even in cases of total edentulism

MOUNTING INSIDE THE MOUTH WITH ANTAGONIST



OPTICAL SCANNING IN LABORATORY



RADIOLOGICAL GUIDE DELIVERY TO THE PATIENT FOR THE CT/CBCT EXAM



- Arch to operate, antagonist and centric (raised)
- Material: silicone or alginate
- Maximum extension, non-functionalised (anatomical)

Reproduce the exact situation of the mouth at the time of surgery

- Extra-hard plaster model casting (no imperfections or bubbles)

The procedure is comparable to that used for the construction of a total prosthesis

The duplicate of the aesthetic test must be precise and with transparent and NOT radiopaque material (orthodontic resin type), after having restored the vestibular flange

COMPATIBLE PROTOCOL

Although it is NOT recommended, if the patient has adequate removable prosthesis (from the aesthetic and functional aspect), proceed with its duplication after having relined it.

For the radiological guide to be stabilised, position the bite in the patient's mouth and reline towards the antagonist with the same silicone supplied with the Evobite kit, taking care that you keep the bite steady

ALTERNATIVE PROTOCOL

Mounting the Evobite to the antagonist can also be performed in the articulator ONLY IF a RAISED centric was used during the initial phase. The same should be done in the case of guided surgery on two arches during the same operation

NOTE: THE COMPONENTS TO CONSTRUCT THE EVOBITE ARE DISPOSABLE