


Sleeves Pricelist

NOTES:

1. 3DIEMME RealGUIDE universal guided surgery kit is used with **3DM00604CAD1.10** sleeve;
2. **3DM00604CADxx** sleeves have RealGUIDE 3DM00604CAD1.10 same outer shape, but they are different in internal diameter (D_{int}) and height. They are dedicated to the original guided surgery kit you can read in each sleeve description;
3. **3DM00605xxx** and **3DM00609xxx** sleeves are comparable to the original sleeves of each manufacturer guided surgery kit.
4. All prices are VAT excluded.


ARTICLE CODE	DESCRIPTION	PRICE [€]
3DM00604CAD0	Surgical guide sleeve RG CAD D_{int} = 5.05 mm Compatible kit: ALLMED	5
3DM00604CAD1.10	Surgical guide sleeve RG CAD D_{int} = 5.05 mm Compatible kit: 3DIEMME RG, ZB Z3D, Nobelguide RP	5
3DM00604CAD1.SMALL	Surgical guide sleeve RG CAD D_{int} = 3.55 mm Compatible kit: 3DIEMME RG MINI	5
3DM00604CAD10	Surgical guide sleeve RG CAD D_{int} = 5.3 mm Compatible kit: DIO Implant	5
3DM00604CAD11	Surgical guide sleeve RG CAD D_{int} = 4.9 mm Compatible kit: Ankylos WD	5
3DM00604CAD2	Surgical guide sleeve RG CAD D_{int} = 4.65 mm Compatible kit: Astra Facilitate MS1	5
3DM00604CAD3	Surgical guide sleeve RG CAD D_{int} = 4.5 mm Compatible kit: Xive ND, Ankylos ND, CAMLOG	5
3DM00604CAD4	Surgical guide sleeve RG CAD D_{int} = 4.2 mm Compatible kit: Anthogyr MS1	5
3DM00604CAD5	Surgical guide sleeve RG CAD D_{int} = 5.2 mm Compatible kit: Xive WD, Anthogyr MS2	5
3DM00604CAD6	Surgical guide sleeve RG CAD D_{int} = 4.1 mm Compatible kit: Nobelguide NP	5
3DM00604CAD8	Surgical guide sleeve RG CAD D_{int} = 5.1 Compatible kit: Navigator 3I LARGE	5
3DM00604CAD9	Surgical guide sleeve RG CAD D_{int} = 4.7 mm	5





ARTICLE CODE	DESCRIPTION	PRICE [€]
3DM006053iLARGE	Surgical guide sleeve - Kit Navigator 3i LARGE	5
3DM006053iSMALL	Surgical guide sleeve - Kit Navigator 3i SMALL	5
3DM00605ALPHASLS	Surgical guide sleeve - Kit AlphaBio SMALL	20
3DM00605ALPHASLL	Surgical guide sleeve - Kit AlphaBio LARGE	20
3DM00605AST.EV.ND	Surgical guide sleeve - Kit Astra EV ND SMALL	5
3DM00605AST.EV.WD	Surgical guide sleeve - Kit Astra EV WD LARGE	5
3DM00605BIOH38	Surgical guide sleeve - Kit BioHorizons 3.8 mm	10
3DM00605BIOH46	Surgical guide sleeve - Kit BioHorizons 4.6 mm	10
3DM00605BIOH58	Surgical guide sleeve - Kit BioHorizons 5.8 mm	10
3DM00605CAM33	Surgical guide sleeve - Kit CAMLOG 3.3 mm	8
3DM00605CAM38	Surgical guide sleeve - Kit CAMLOG 3.8 mm	8
3DM00605CAM43	Surgical guide sleeve - Kit CAMLOG 4.3 mm	8
3DM00605F4GS	Surgical guide sleeve - Kit Fixo Biomec F4GS	5
3DM00605FMDSM	Surgical guide sleeve - Kit FMD Small	5
3DM00605FMDLG	Surgical guide sleeve - Kit FMD Large	5
3DM00605GEASS42	Surgical guide sleeve - Kit Geass 4.2 mm	10
3DM00605GEASS52	Surgical guide sleeve - Kit Geass 5.2 mm	10
3DM00605GHISML	Surgical guide sleeve - Kit Ghimas 4.15 mm (small)	5
3DM00605GHILRG	Surgical guide sleeve - Kit Ghimas 5.50 mm (large)	5
3DM00605INTLK	Surgical guide sleeve - Kit INTRA-LOCK (PEEK)	20
3DM00605ISOSM	Surgical guide sleeve - Kit ISOMED Small	5
3DM00605ISOLG	Surgical guide sleeve - Kit ISOMED Large	5
3DM00605iRES	Surgical guide sleeve - Kit iRES	5
3DM00605MGN	Surgical guide sleeve - Kit Megagen	5
3DM00605MGNWD	Surgical guide sleeve - Kit Megagen Wide	5
3DM00605MIS	Surgical guide sleeve - Kit MGUIDE MIS H4	5
3DM00605MISNRW	Surgical guide sleeve - Kit MGUIDE MIS Narrow H4	5
3DM00605NBL.NP	Surgical guide sleeve - Kit Nobel NP	5
3DM00605NBL.RP	Surgical guide sleeve - Kit Nobel RP	5
3DM00605NBL.WP	Surgical guide sleeve - Kit Nobel WP	5
3DM00605NEOSS	Surgical guide sleeve - Kit Neoss	5
3DM00605SM415.EX	Surgical guide sleeve - Kit Sweden & Martina GS-415 EX	5
3DM00605SM550.EX	Surgical guide sleeve - Kit Sweden & Martina GS-550 EX	5
3DM00605STR	Surgical guide sleeve - Kit Straumann	5





3DM00608CAD	Surgical guide pin sleeve - Kit 3DIEMME RealGUIDE	5
3DM00609ALPHASLSE	Surgical guide pin sleeve - Kit AlphaBio	20
3DM00609FMD	Surgical guide pin sleeve - Kit FMD	5
3DM00609GEASS	Surgical guide pin sleeve - Kit GEASS	15
3DM00609ISO	Surgical guide pin sleeve - Kit ISOMED	5
3DM00609MIS	Surgical guide pin sleeve - Kit MGUIDE MIS	5

Rev.01 - 2020



3DIEMME S.r.l.
P.IVA 02938110133
www.3diemme.it
support@3diemme.it

Sede Legale:
Via Risorgimento,9
22063 Cantù – (CO)
Italia

Sede Operativa:
Via Delle Baraggette,10
22060 Figino Serenza – (CO)
Italia

Telefono:
+39 031 7073353